



INVERTED MICROSCOPES

Model ZEN-INVI

For tissue culture examination meant for Embryology, Immunology, Pharmacology, Microbiology, Virology Water management, Biochemistry, Biotechnology and more...

SPECIFICATIONS

Model	ZEN-INVI
Head	Binocular providing adequate comfort viewing, diopter adjustment ring on one tube are standard feature to facilitate compensation. Graduated interpupillary distance adjustments is from 52mm to 76mm.
Objectives	LWD plan 4x/16mm, 10x/11.2mm, 20x/6.02mm.
Eyepieces	Wide Field 10x paired Field number 18mm
Stage	Fixed plain stage to accomodate large size culture disc.
Focusing	Coarse and fine focussing mechanism is designed for frequent use provides a smooth focus adjustments.
Condenser	N.A 0.3 long working distance with iris diaphragm and filter holder adjustable to make enough space above the stage.
Illumination	Built-in 6v-20w halogen lamp with intensity regulator, input 220-230v AC 50/60Hz.
Standard accessories	One spare lamp, 2 rubber eye guard Blue filter, instruction manual.
Packing	Microscope packed in plywood case with lock & key.