



## STEREOSCOPIC BINOCULAR MICROSCOPE



# OLYMPUS®

		<b>MS13</b>	<b>MS24</b>
• <b>Body :</b>	<b>Objectives</b>	1X & 3X	2X & 4X
	<b>Working distance</b>	105mm	
	<b>Magnification</b>	10x / 30x (10x eyepiece) 15x / 45x (15x eyepiece)	20x / 40x (10x eyepiece) 30x / 60x (15x eyepiece)
• <b>Eyepieces :</b>	<b>High eyepoint Super widefield eyepieces SWH10X with Field Number 23</b>		
• <b>Field of view :</b>	<b>Objective</b>	1X    3X	2X    4X
	<b>Visual field (mm)</b>	23    7.7	11.5    5.7
• <b>Light source :</b>	<b>Reflected</b>	6V 10W lamp with adjustable brightness	
	<b>Transmitted</b>	5W fluorescent lamp	

## ZENITH ENGINEERS

No. 4/5 Ram Babu Verma Marg, Civil Lines,  
Agra - 282002, India.

[www.BioZen.co.in](http://www.BioZen.co.in)

ph: (0562) 2523036

fx: (0562) 2852141

em: [mohit@BioZen.co.in](mailto:mohit@BioZen.co.in)